LAB2-AC2 Group E - Calubaquib, Reyes, Sanguyo

## [Francis Sanguyo's copy]

# **Objectives:**

This experiment aims to:

- 1. extract the essential oils (EOs) present in fresh lemongrass via steam distillation, and
- 2. determine the mass spectra, possible chemical structures, and relative abundance of major components of lemongrass EOs using GC-MS with electron ionization (EI).


### Reagents

Reagent Name and Structure	Physical and Chemical Properties	Safety Precautions	
Distilled Water (H <sub>2</sub> O)  H O	- colorless, odorless liquid - MW = 18.02 g/mol	Non-irritant, generally safe to handle  Disposal: sink	
Lemongrass Sample	- solid sample - contains liquid oil which may be extracted via steam distillation	- clean lemongrass samples before and after handling - handle lemongrass samples with washed, clean hands	
Anhydrous Sodium Sulfate (Na <sub>2</sub> SO <sub>4</sub> ) $2Na^{+} \begin{bmatrix} O \\ I \\ O \end{bmatrix}^{2-}$	<ul> <li>white solid with characteristic odor</li> <li>hygroscopic; absorbs surrounding moisture</li> <li>clumps together in solution after absorbing moisture</li> <li>soluble in water</li> <li>MW = 142.04 g/mol</li> </ul>	- chemical is not considered hazardous by 2012 OSHA  Disposal: <u>D499</u>	
n-Hexane (C <sub>6</sub> H <sub>14</sub> )	- clear, colorless liquid with characteristic odor - non-polar solvent; dissolves many oils - MW = 86.18 g/mol - volatile and low BP (69 °C); serves as solvent for GC	- causes skin irritation - may be fatal if it is swallowed or if it enters airways - toxic to aquatic life with long-lasting effects  Disposal: G704	

Acetone (C <sub>3</sub> H <sub>6</sub> O) O CH <sub>3</sub>	<ul> <li>clear, colorless liquid with pungent odor resembling of nail polish</li> <li>polar aprotic solvent; dissolves many oils</li> <li>MW = 58.08 g/mol</li> <li>volatile and low BP (56 °C); serves as solvent for GC</li> </ul>	<ul> <li>causes serious eye irritation</li> <li>may cause drowsiness or dizziness</li> <li>highly flammable liquid and vapor</li> <li>Disposal: <u>G704</u></li> </ul>
Helium Gas (He)	- colorless, odorless gas - MW = 4 g/mol - serves as mobile phase for GC	- may displace oxygen and cause suffocation - contains gas under pressure; may explode if heated

### **Equipment:**

- Glassware for Steam Distillation Setup
- Water Circulating Pump
- Heating Mantle
- Pasteur pipet
- Gas Cylinder for Helium Carrier Gas
- GC-MS Setup



#### **Procedure:**

- I. Extraction of EOs from Lemongrass via Steam Distillation
  - a. Weigh 50 g of fresh lemongrass. Carefully cut the sample into 1-cm-long bits.
  - b. Transfer the lemongrass sample to a 1 L round bottom flask (RBF).
  - c. Add 500 mL distilled water
  - d. Assemble the steam distillation setup, as shown in the figure above.
  - e. Turn on the water circulating pump. Make sure there are no leaks.
  - f. Place the RBF on the heating mantle, and turn on the heating.
  - g. <u>Distill the lemongrass mixture for two hours.</u>
  - h. Observe carefully the appearance of the distillate (aqueous layer + top oil layer)

- i. Once finished with distillation, turn off the heating mantle.
- j. Collect lemongrass EO from the top oil layer by using a Pasteur pipette and sample vial.
- k. Add anhydrous Na<sub>2</sub>SO<sub>4</sub> to the EO samples and allow the oil to dry.
- 1. Store the oil at a cold temperature of 4 °C before GC-MS analysis.

### II. GC-MS Analysis of Lemongrass EOs

- a. Estimate the amounts of solvent required to dilute lemongrass EO by 1000-fold.
- b. Dilute the EO by 1000-fold using a selected solvent (may be hexane, acetone, etc.)
- c. Turn on the following in order: Helium carrier gas, GC system, MS system, and finally the assigned computer for GC-MS.
- d. Perform auto start-up to initialize the vacuum pump and ion source temperature.
- e. <u>Leave the vacuum to stabilize for at least two hours.</u>
- f. Perform leak change and functioning.
- g. Open a GC-MS method, and download the method.
- h. Inject solvent headspace to determine the solvent cut-off time.
- i. Analyze the essential oil solution that was prepared.

# Literature Data

I. Notable Components Found in Lemongrass EOs

Table 1. Components of Lemongrass Oil Extract According to Carlson et al. (2001)

ruste 1. Comp	Thenes of Lemongrass on I	Extract According to Caris	311 ct at. (2001)
Component Name	Structure	Area% at 90 bar and 23 °C	Area% at 85 bar and 50 °C
Gerenial	сно	49.79	51.07
Neral	сно	27.47	31.32
Myrcene	CH <sub>2</sub> CH <sub>2</sub> CH <sub>2</sub>	5.29	5.01
Geraniol	CH <sub>3</sub> OH H <sub>3</sub> C CH <sub>3</sub>	2.31	1.82
Pentacosane	(CH <sub>2</sub> ) <sub>23</sub> CH <sub>3</sub>	1.06	0.13
Linalool	H <sub>3</sub> C OH CH <sub>2</sub> H <sub>3</sub> C CH <sub>3</sub>	0.96	1.27
Heptacosane	(CH <sub>2</sub> ) <sub>25</sub> CH <sub>3</sub> CH <sub>3</sub>	0.75	0.16
Limonene		trace	trace

*Table 2.* Components of Lemongrass Essential Oil in Underground and Above-Ground Shoots According to Kieltyka-Dadasiewicz *et al.* (2021)

Component Name	Structure	Relative Percentage in Underground Shoots	Relative Percentage in Above-Ground Shoots
Gerenial	сно	41.1	48.0
Neral	сно	31.3	34.6
Elemol	HO HO	14.0	0.8
Geranyl Acetate	CH <sub>3</sub> O CH <sub>3</sub>	2.1	0.5
6-methyl-5-heptene -2-one	CH <sub>3</sub> O CH <sub>3</sub>	1.3	1.6
Geraniol	CH <sub>3</sub> OH H <sub>3</sub> C CH <sub>3</sub>	1.2	0.4
Linalool	H <sub>3</sub> C OH CH <sub>2</sub> CH <sub>3</sub>	0.5	0.7
Limonene		0.2	0.4

# II. Reference Tables for MS Fragmentation

Table 3. Summary of Detected Fragments and Lost Fragments in MS Spectra of Organic Compounds

Decrease in m/z by	Radical lost	Class of compound
1	H-	any
14	NH <sub>2</sub> ·	amine
15	CH3-	methyl
17	OH-	alcohol
28	co.	carbonyl
29	CH3CH2-	ethyl
29	сно-	aldehyde
31	CH₂OH-	primary alcohol
31	OCH3-	ester
43	COCH₃·	ketone
43	CH3CH2CH2-	propyl
45	соон-	carboxyl
77	C6H5.	phenyl

Table 4. Summary of Mechanisms Involved in Fragment Production (Dunnivant, 2017; McMurry, 2016)

#	Functional Group	Mechanism	Example
1	Branched Alkane	Sigma Cleavage	$\begin{bmatrix} CH_3 \\ H_3C - C - CH_3 \\ CH_3 \\ CH_3 \end{bmatrix}^{+} \cdot \begin{bmatrix} CH_3 \\ H_3C - C^+ \\ CH_3 \\ CH_3 \end{bmatrix}$
2	Alkene	Allylic Cleavage	CH <sub>CH</sub> CH <sub>2</sub>
3	Alcohol	Alpha Cleavage	$\begin{bmatrix} RCH_2 & OH \\ C & C \end{bmatrix}^{+} \xrightarrow{Cleavage}  RCH_2 \cdot  +  \begin{bmatrix} \vdots \ddot{OH} & & \dagger \ddot{OH} \\ C + & C & & \\ C & C \end{bmatrix}$
4	Alcohol	Dehydration	$\begin{bmatrix} H & OH \\ C - C \end{bmatrix}^{+} \xrightarrow{\text{Dehydration}} H_2O + \begin{bmatrix} C = C \end{bmatrix}^{+}$
5	Aldehyde or Ketone	Alpha Cleavage	$\begin{bmatrix} \begin{matrix} O \\ \parallel \\ C \\ R \end{matrix} \end{bmatrix}^{+ \cdot} \xrightarrow{\text{Alpha}}  R \cdot  +  \begin{bmatrix} :O: & :O^{+} \\ \parallel \\ C^{+} & \longleftrightarrow & C \\ \mid \\ R' & & R' \end{bmatrix}$
6	Aldehyde or Ketone	McLafferty Rearrangement	McLafferty rearrangement C +   C C   C   C   C   C   C   C   C   C
7	Carboxyl Ester	Alpha Cleavage	$ \begin{array}{cccccccccccccccccccccccccccccccccccc$
8	Carboxyl Ester	McLafferty Rearrangement	H to: H <sub>2</sub> C C

**Data and Observations:** 

#### **References:**

- Carlson, L., Machado, R., Spricago, C., Pereira, L., and Bolzan, A. (2001). Extraction of lemongrass essential oil with dense carbon dioxide. *Journal of Supercritical Fluids*, *21*, 33-39. https://doi.org/10.1016/S0896-8446(01)00085-7
- Dunnivant, F. (2017). *Fragmentation and Interpretation of Spectra*. A basic introduction to Gas Chromatography, Liquid Chromatography, Capillary Electrophoresis, and Mass Spectrometry. https://sites.google.com/whitman.edu/gclccems/home
- Kieltya-Dadasiewicz, A., Ludwiczuk, A., Teraseviciene, Z., Michalak, M., Glowacka, A., Baj, T., Krecisz, and Krochmal-Marczak, B. (2021). Chemical and Nutritional Compounds of Different Parts of Lemongrass (Cymbopogon citratus (DC) Stapf.) Cultivated in Temperate Climate of Poland. *Journal of Oleo Science*, 70(1), 125-133.

McMurry, J. (2016). Organic Chemistry (9th ed.). Cengage Learning.

 $\frac{\text{fishersci.com/store/msds?partNumber=S421500\&productDescription=SODIUM+SULFATE+CERT+AC}{\text{S+500G\&vendorId=VN00033897\&countryCode=US\&language=en}}$ 

chemos.de/import/data/msds/GB en/110-54-3-A0001755-GB-en.pdf

sigmaaldrich.com/PH/en/sds/sigald/179124?userType=undefined

airgas.com/msds/001025.pdf